

Ubiquitin-protein ligase E3A

Catalog No: #AP88685

Package Size: #AP88685-1 50ug #AP88685-2 100ug #AP88685-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ubiquitin-protein ligase E3A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	E6AP ubiquitin-protein ligase,HECT-type ubiquitin transferase E3A,Human papillomavirus E6-associated protein,Oncogenic protein-associated protein E6-AP,Renal carcinoma antigen NY-REN-54,E6AP,EPVE6AP,HPVE6A
Accession No.	Uniprot ID:Q05086
Uniprot	Q05086
GeneID	7337;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:UBE3A

Note: This product is for in vitro research use only