

Non-receptor tyrosine-protein kinase TNK1

Catalog No: #AP88688



Package Size: #AP88688-1 50ug #AP88688-2 100ug #AP88688-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Non-receptor tyrosine-protein kinase TNK1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Kinase of embryonic stem cells,Kos1
Accession No.	Uniprot ID:Q99ML2
Uniprot	Q99ML2
GeneID	83813;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TNK1

Note: This product is for in vitro research use only