

TATA box-binding protein-like protein 1

Catalog No: #AP88690



Package Size: #AP88690-1 50ug #AP88690-2 100ug #AP88690-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	TATA box-binding protein-like protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	TBP-like protein 1, TATA box-binding protein-related factor 2, TBP-related factor 2, TBP-like factor, TBP-related protein, TLF, TLP, TRF2, TRP
Accession No.	Uniprot ID:Q32LB1
Uniprot	Q32LB1
GeneID	613818;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0) gene_full_name: TBPL1

Note: This product is for in vitro research use only