

Pro-thyrotropin-releasing hormone

Catalog No: #AP88691



Package Size: #AP88691-1 50ug #AP88691-2 100ug #AP88691-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pro-thyrotropin-releasing hormone
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Pro-TRH,Prothyroliberin
Accession No.	Uniprot ID:P01150
Uniprot	P01150
GeneID	25569;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Trh

Note: This product is for in vitro research use only