

UDP-glucuronosyltransferase 1-6

Catalog No: #AP88693



Package Size: #AP88693-1 50ug #AP88693-2 100ug #AP88693-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	UDP-glucuronosyltransferase 1-6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	UDPGT 1-6, Phenol UDP-glucuronosyltransferase, UDP-glucuronosyltransferase 1A6, UGT1A6, UGP1A1, UGT1A7, Ugt1, Ugt1a6a, Ugt1a7
Accession No.	Uniprot ID:Q64435
Uniprot	Q64435
GeneID	94284;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ugt1a6

Note: This product is for in vitro research use only