

Thrombomodulin

Catalog No: #AP88695



Package Size: #AP88695-1 50ug #AP88695-2 100ug #AP88695-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Thrombomodulin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	TM,Fetomodulin,CD141
Accession No.	Uniprot ID:P06579
Uniprot	P06579
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:THBD

Note: This product is for in vitro research use only