

# TGF-beta-activated kinase 1 and MAP3K7-binding protein 2

Catalog No: #AP88696

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP88696-1 50ug #AP88696-2 100ug #AP88696-3 1mg

## Description

Product Name	TGF-beta-activated kinase 1 and MAP3K7-binding protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Mitogen-activated protein kinase kinase kinase 7-interacting protein 2,TGF-beta-activated kinase 1-binding protein 2,Map3k7ip2
Accession No.	Uniprot ID:Q5U303
Uniprot	Q5U303
GeneID	308267;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Tab2

Note: This product is for in vitro research use only