

Thyroxine-binding globulin

Catalog No: #AP88699



Package Size: #AP88699-1 50ug #AP88699-2 100ug #AP88699-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Thyroxine-binding globulin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Serpin A7,T4-binding globulin,TBG
Accession No.	Uniprot ID:Q9TT36
Uniprot	Q9TT36
GeneID	282518;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SERPINA7

Note: This product is for in vitro research use only