

Signal recognition particle 19 kDa protein

Catalog No: #AP88700



Package Size: #AP88700-1 50ug #AP88700-2 100ug #AP88700-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Signal recognition particle 19 kDa protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	SRP19
Accession No.	Uniprot ID:P09132
Uniprot	P09132
GeneID	6728;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SRP19

Note: This product is for in vitro research use only