

SteryI-sulfatase

Catalog No: #AP88704



Package Size: #AP88704-1 50ug #AP88704-2 100ug #AP88704-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SteryI-sulfatase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Arylsulfatase C,ASC,Steroid sulfatase,SteryI-sulfate sulfohydrolase,ARSC1
Accession No.	Uniprot ID:P08842
Uniprot	P08842
GeneID	412;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:STS

Note: This product is for in vitro research use only