

V-type proton ATPase catalytic subunit A

Catalog No: #AP88706



Package Size: #AP88706-1 50ug #AP88706-2 100ug #AP88706-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	V-type proton ATPase catalytic subunit A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	V-ATPase 116 kDa isoform a2,Vacuolar proton translocating ATPase 116 kDa subunit a isoform 2,ATP6N1B
Accession No.	Uniprot ID:O97681
Uniprot	O97681
GeneID	338038;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ATP6V1A

Note: This product is for in vitro research use only