

## Vascular endothelial growth factor receptor 3

Catalog No: #AP88708



Package Size: #AP88708-1 50ug #AP88708-2 100ug #AP88708-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Vascular endothelial growth factor receptor 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	VEGFR-3,Fms-like tyrosine kinase 4,FLT-4,Tyrosine-protein kinase receptor FLT4,VEGFR3
Accession No.	Uniprot ID:P35916
Uniprot	P35916
GeneID	2324;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:FLT4

Note: This product is for in vitro research use only