

## Vasoactive intestinal polypeptide receptor 2

Catalog No: #AP88713



Package Size: #AP88713-1 50ug #AP88713-2 100ug #AP88713-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Vasoactive intestinal polypeptide receptor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	VIP-R-2, Helodermin-preferring VIP receptor, Pituitary adenylate cyclase-activating polypeptide type III receptor, PACAP type III receptor, VPAC2, VIP2R
Accession No.	Uniprot ID:P41587
Uniprot	P41587
GeneID	7434;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: VIPR2

Note: This product is for in vitro research use only