

# Visinin-like protein 1

Catalog No: #AP88714



Package Size: #AP88714-1 50ug #AP88714-2 100ug #AP88714-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Visinin-like protein 1                                      |
| Brief Description  | Recombinant Protein   |
| Host Species       | E.coli  |
| Purification       | Greater than 90% by SDS-PAGE                                |
| Species Reactivity | Rt  |
| Other Names        | VILIP,21 kDa CABP,Neural visinin-like protein 1,NVL-1,Visl1 |
| Accession No.      | Uniprot ID:P62762   |
| Uniprot            | P62762  |
| Tag Info           | His   |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)      |

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Vsnl1

Note: This product is for in vitro research use only