

Wiskott-Aldrich syndrome protein

Catalog No: #AP88717



Package Size: #AP88717-1 50ug #AP88717-2 100ug #AP88717-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Wiskott-Aldrich syndrome protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	WASp,IMD2
Accession No.	Uniprot ID:P42768
Uniprot	P42768
GeneID	7454;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:WAS

Note: This product is for in vitro research use only