

## Erythrocyte band 7 integral membrane protein

Catalog No: #AP88723



Package Size: #AP88723-1 50ug #AP88723-2 100ug #AP88723-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Erythrocyte band 7 integral membrane protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Protein 7.2b,Stomatin,Epb7.2,Epb72
Accession No.	Uniprot ID:P54116
Uniprot	P54116
GeneID	13830;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Stom

Note: This product is for in vitro research use only