

Transmembrane protease serine 4

Catalog No: #AP88724



Package Size: #AP88724-1 50ug #AP88724-2 100ug #AP88724-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transmembrane protease serine 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Channel-activating protease 2,mCAP2,Cap2
Accession No.	Uniprot ID:Q8VCA5
Uniprot	Q8VCA5
GeneID	214523;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tmprss4

Note: This product is for in vitro research use only