

Tumor necrosis factor receptor superfamily member 10B

Catalog No: #AP88725

Package Size: #AP88725-1 50ug #AP88725-2 100ug #AP88725-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor receptor superfamily member 10B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Death receptor 5,MK,CD262,Dr5,Killer
Accession No.	Uniprot ID:Q9QZM4
Uniprot	Q9QZM4
GeneID	21933;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tnfrsf10b

Note: This product is for in vitro research use only