

tRNA (adenine(58)-N(1))-methyltransferase, mitochondrial

Catalog No: #AP88727

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP88727-1 50ug #AP88727-2 100ug #AP88727-3 1mg

Description

Product Name	tRNA (adenine(58)-N(1))-methyltransferase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Accession No.	Uniprot ID:Q9BVS5
Uniprot	Q9BVS5
GeneID	55006;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: TRMT61B

Note: This product is for in vitro research use only