

## Vacuole membrane protein 1

Catalog No: #AP88735



Package Size: #AP88735-1 50ug #AP88735-2 100ug #AP88735-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Vacuole membrane protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Transmembrane protein 49,TDC1,TMEM49,HSPC292
Accession No.	Uniprot ID:Q96GC9
Uniprot	Q96GC9
GeneID	81671;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:VMP1

Note: This product is for in vitro research use only