

Transcription factor E2-alpha

Catalog No: #AP88736



Package Size: #AP88736-1 50ug #AP88736-2 100ug #AP88736-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transcription factor E2-alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Class B basic helix-loop-helix protein 21,bHLHb21,Immunoglobulin enhancer-binding factor E12/E47,Immunoglobulin transcription factor 1,Kappa-E2-binding factor,Transcription factor 3,TCF-3,Transcription factor ITF-1,BHLHB21,E2A,ITF1
Accession No.	Uniprot ID:P15923
Uniprot	P15923
GeneID	6929;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TCF3

Note: This product is for in vitro research use only