

Mitochondrial uncoupling protein 2

Catalog No: #AP88749



Package Size: #AP88749-1 50ug #AP88749-2 100ug #AP88749-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mitochondrial uncoupling protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	UCP 2,Solute carrier family 25 member 8,SLC25A8
Accession No.	Uniprot ID:Q3SZI5
Uniprot	Q3SZI5
GeneID	281562;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:UCP2

Note: This product is for in vitro research use only