

Tryptophan--tRNA ligase, cytoplasmic

Catalog No: #AP88753



Package Size: #AP88753-1 50ug #AP88753-2 100ug #AP88753-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tryptophan--tRNA ligase, cytoplasmic
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Tryptophanyl-tRNA synthetase, TrpRS, Wrs
Accession No.	Uniprot ID:P32921
Uniprot	P32921
GeneID	22375;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Wrs

Note: This product is for in vitro research use only