

## Tsukushin

Catalog No: #AP88754



Package Size: #AP88754-1 50ug #AP88754-2 100ug #AP88754-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tsukushin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Tsukushi, Early insulin-induced hepatic gene protein, EIIH, Leucine-rich repeat-containing protein 54, Eiih, Lrrc54, Tsk
Accession No.	Uniprot ID: Q6QMY6
Uniprot	Q6QMY6
GeneID	308843;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Tsku

Note: This product is for in vitro research use only