

Triosephosphate isomerase

Catalog No: #AP88756



Package Size: #AP88756-1 50ug #AP88756-2 100ug #AP88756-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Triosephosphate isomerase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	TIM, Triose-phosphate isomerase, Tpi
Accession No.	Uniprot ID:P17751
Uniprot	P17751
GeneID	21991;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Tpi1

Note: This product is for in vitro research use only