

# RISC-loading complex subunit TARBP2

Catalog No: #AP88757



Package Size: #AP88757-1 50ug #AP88757-2 100ug #AP88757-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	RISC-loading complex subunit TARBP2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Accession No.	Uniprot ID:Q3SWU0
Uniprot	Q3SWU0
GeneID	363006;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Tarbp2

Note: This product is for in vitro research use only