

## E3 ubiquitin-protein ligase TRIM22

Catalog No: #AP88758



Package Size: #AP88758-1 50ug #AP88758-2 100ug #AP88758-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	E3 ubiquitin-protein ligase TRIM22
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	50 kDa-stimulated trans-acting factor,RING finger protein 94,RING-type E3 ubiquitin transferase TRIM22,Staf-50,Tripartite motif-containing protein 22,RNF94,STAF50
Accession No.	Uniprot ID:Q8IYM9
Uniprot	Q8IYM9
GeneID	10346;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:TRIM22

Note: This product is for in vitro research use only