

Lymphotoxin-alpha

Catalog No: #AP88759



Package Size: #AP88759-1 50ug #AP88759-2 100ug #AP88759-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Lymphotoxin-alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	LT-alpha,TNF-beta,Tumor necrosis factor ligand superfamily member 1,TNFB,TNFSF1
Accession No.	Uniprot ID:Q06600
Uniprot	Q06600
GeneID	280845;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LTA

Note: This product is for in vitro research use only