

Tumor necrosis factor receptor superfamily member 13B



Catalog No: #AP88760

Orders: order@signalwayantibody.com

Package Size: #AP88760-1 50ug #AP88760-2 100ug #AP88760-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor receptor superfamily member 13B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Transmembrane activator and CAML interactor,CD267,TACI
Accession No.	Uniprot ID:O14836
Uniprot	O14836
GeneID	23495;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TNFRSF13B

Note: This product is for in vitro research use only