

Vang-like protein 2

Catalog No: #AP88763



Package Size: #AP88763-1 50ug #AP88763-2 100ug #AP88763-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Vang-like protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Loop-tail protein 1 homolog, Strabismus 1, Van Gogh-like protein 2, KIAA1215, STB1
Accession No.	Uniprot ID:Q9ULK5
Uniprot	Q9ULK5
GeneID	57216;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: VANGL2

Note: This product is for in vitro research use only