

Signal transducer and activator of transcription 5A

Catalog No: #AP88768



Package Size: #AP88768-1 50ug #AP88768-2 100ug #AP88768-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Signal transducer and activator of transcription 5A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Mammary gland factor,Mgf,Mpf
Accession No.	Uniprot ID:P42230
Uniprot	P42230
GeneID	20850;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Stat5a

Note: This product is for in vitro research use only