

Tumor necrosis factor ligand superfamily member 13B

Catalog No: #AP88773



Package Size: #AP88773-1 50ug #AP88773-2 100ug #AP88773-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor ligand superfamily member 13B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	B-cell-activating factor,BAFF,CD257,Baff
Accession No.	Uniprot ID:Q9WU72
Uniprot	Q9WU72
GeneID	24099;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Tnfsf13b

Note: This product is for in vitro research use only