

## Tapasin-related protein

Catalog No: #AP88777



Package Size: #AP88777-1 50ug #AP88777-2 100ug #AP88777-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tapasin-related protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	TAPASIN-R, TAP-binding protein-like, TAP-binding protein-related protein, TAPBP-R, Tapasin-like
Accession No.	Uniprot ID:Q8VD31
Uniprot	Q8VD31
GeneID	213233;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Tapbp1

Note: This product is for in vitro research use only