

Protein S100-A6

Catalog No: #AP88780



Package Size: #AP88780-1 50ug #AP88780-2 100ug #AP88780-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein S100-A6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Calcyclin,S100 calcium-binding protein A6
Accession No.	Uniprot ID:Q98953
Uniprot	Q98953
GeneID	373951;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:S100A6

Note: This product is for in vitro research use only