

# Triosephosphate isomerase

Catalog No: #AP88784



Package Size: #AP88784-1 50ug #AP88784-2 100ug #AP88784-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Triosephosphate isomerase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	TIM, Triose-phosphate isomerase, TPI
Accession No.	Uniprot ID:P60174
Uniprot	P60174
GeneID	7167;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: TPI1

Note: This product is for in vitro research use only