

## Trimethyllysine dioxygenase, mitochondrial

Catalog No: #AP88785



Package Size: #AP88785-1 50ug #AP88785-2 100ug #AP88785-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Trimethyllysine dioxygenase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Epsilon-trimethyllysine 2-oxoglutarate dioxygenase,TML hydroxylase,TML-alpha-ketoglutarate dioxygenase,TML dioxygenase,RCJMB04_1d18
Accession No.	Uniprot ID:Q5F4B3
Uniprot	Q5F4B3
GeneID	422296;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:TMLHE

Note: This product is for in vitro research use only