

mRNA decay activator protein ZFP36

Catalog No: #AP88788



Package Size: #AP88788-1 50ug #AP88788-2 100ug #AP88788-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	mRNA decay activator protein ZFP36
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Tristetraprolin, TTP, Zinc finger protein 36, Zfp-36
Accession No.	Uniprot ID:P53781
Uniprot	P53781
GeneID	282127;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: ZFP36

Note: This product is for in vitro research use only