

Proto-oncogene Wnt-1

Catalog No: #AP88790



Package Size: #AP88790-1 50ug #AP88790-2 100ug #AP88790-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Proto-oncogene Wnt-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Proto-oncogene Int-1,Int-1,Wnt-1
Accession No.	Uniprot ID:P04426
Uniprot	P04426
GeneID	22408;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Wnt1

Note: This product is for in vitro research use only