

Protein Wnt-11

Catalog No: #AP88791



Package Size: #AP88791-1 50ug #AP88791-2 100ug #AP88791-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein Wnt-11
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	WNT-11
Accession No.	Uniprot ID:P49339
Uniprot	P49339
GeneID	395562;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:WNT11

Note: This product is for in vitro research use only