

## WNT1-inducible-signaling pathway protein 2

Catalog No: #AP88795



Package Size: #AP88795-1 50ug #AP88795-2 100ug #AP88795-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	WNT1-inducible-signaling pathway protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	WISP-2,CCN family member 5,Connective tissue growth factor-like protein,CTGF-L,Ccn5,Ctgfl
Accession No.	Uniprot ID:Q9Z0G4
Uniprot	Q9Z0G4
GeneID	22403;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Wisp2

Note: This product is for in vitro research use only