

Tumor necrosis factor ligand superfamily member 6

Catalog No: #AP88798



Package Size: #AP88798-1 50ug #AP88798-2 100ug #AP88798-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor ligand superfamily member 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CD95 ligand,CD95-L,Fas antigen ligand,Fas ligand,CD178,Apt1lg1,Cd95l,Fasl,gld,Tnfsf6
Accession No.	Uniprot ID:P41047
Uniprot	P41047
GeneID	14103;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Faslg

Note: This product is for in vitro research use only