

Transient receptor potential cation channel subfamily V member 1

Catalog No: #AP88799

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP88799-1 50ug #AP88799-2 100ug #AP88799-3 1mg

Description

Product Name	Transient receptor potential cation channel subfamily V member 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	TrpV1,Osm-9-like TRP channel 1,OTRPC1,Vanilloid receptor 1
Accession No.	Uniprot ID:Q697L1
Uniprot	Q697L1
GeneID	445457;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TRPV1

Note: This product is for in vitro research use only