

Tumor necrosis factor ligand superfamily member 10

Catalog No: #AP88801



Package Size: #AP88801-1 50ug #AP88801-2 100ug #AP88801-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor ligand superfamily member 10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Apo-2 ligand,Apo-2L,TNF-related apoptosis-inducing ligand,Protein TRAIL,CD253,APO2L,TRAIL
Accession No.	Uniprot ID:P50591
Uniprot	P50591
GeneID	8743;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TNFSF10

Note: This product is for in vitro research use only