

Uromodulin

Catalog No: #AP88805



Package Size: #AP88805-1 50ug #AP88805-2 100ug #AP88805-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Uromodulin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Tamm-Horsfall urinary glycoprotein, THP
Accession No.	Uniprot ID:Q91X17
Uniprot	Q91X17
GeneID	22242;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Umod

Note: This product is for in vitro research use only