

V-type proton ATPase catalytic subunit A

Catalog No: #AP88807



Package Size: #AP88807-1 50ug #AP88807-2 100ug #AP88807-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | V-type proton ATPase catalytic subunit A |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Pig |
| Other Names | V-ATPase subunit A,V-ATPase 69 kDa subunit,Vacuolar proton pump subunit alpha,ATP6A1,ATP6V1A1 |
| Accession No. | Uniprot ID:Q29048 |
| Uniprot | Q29048 |
| GeneID | 445531; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ATP6V1A

Note: This product is for in vitro research use only