

# Tryptase beta-2

Catalog No: #AP88809



Package Size: #AP88809-1 50ug #AP88809-2 100ug #AP88809-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Tryptase beta-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Tryptase-2,Mast cell protease 6,rMCP-6,Mcp6,Mcpt6
Accession No.	Uniprot ID:P50343
Uniprot	P50343
GeneID	29268;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Tpsb2

Note: This product is for in vitro research use only