

Tumor necrosis factor receptor superfamily member 6B

Catalog No: #AP88810

Orders: order@signalwayantibody.com

Package Size: #AP88810-1 50ug #AP88810-2 100ug #AP88810-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor receptor superfamily member 6B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Decoy receptor 3,DcR3,Decoy receptor for Fas ligand,M68,DCR3,TR6,UNQ186/PRO212
Accession No.	Uniprot ID:O95407
Uniprot	O95407
GeneID	8771;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TNFRSF6B

Note: This product is for in vitro research use only