

Cystine/glutamate transporter

Catalog No: #AP88812



Package Size: #AP88812-1 50ug #AP88812-2 100ug #AP88812-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cystine/glutamate transporter
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Amino acid transport system xc-, Calcium channel blocker resistance protein CCBR1, Solute carrier family 7 member 11, xCT
Accession No.	Uniprot ID:Q9UPY5
Uniprot	Q9UPY5
GeneID	23657;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: SLC7A11

Note: This product is for in vitro research use only