

Tumor necrosis factor ligand superfamily member 4

Catalog No: #AP88827



Package Size: #AP88827-1 50ug #AP88827-2 100ug #AP88827-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor ligand superfamily member 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	OX40 ligand, OX40L, CD252, Ox40I, Txgp1
Accession No.	Uniprot ID:Q9Z2P3
Uniprot	Q9Z2P3
GeneID	89814;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Tnfsf4

Note: This product is for in vitro research use only