

# Prosalusin

Catalog No: #AP88828



Package Size: #AP88828-1 50ug #AP88828-2 100ug #AP88828-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Prosalusin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Torsin family 2 member A,Torsin-2A
Accession No.	Uniprot ID:P0C7W1
Uniprot	P0C7W1
GeneID	534311;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:TOR2A

Note: This product is for in vitro research use only